


PERFECT WALKER II UNBOX & SET-UP GUIDE

Setting up your Perfect Walker is very simple - just follow these easy instructions. You will be using your walker in just a few minutes. We recommend that you lay out the components to ensure that you have easy access to them at each step.

Page 1

Step 1

Unpack the Perfect Walker II and lay out the parts for easy access.



Step 2

Loosen and remove the 4 large thumbwheel screws (see image below) from the black lower frame.



Step 3

Fold out each of the 4 wheel/leg combinations to their full positions, ensuring that the screw holes are lined up between the frame and the legs.

Step 4

Insert the 4 thumbwheel screws that were removed in Step 2 into the 4 holes in the lower legs. The thumbwheel screws are to be inserted from the outside of the frame. The bolts may be a little tight when inserting into the holes - this is normal. Once inserted as far as they can go, you will screw them into place and lock the wheels into place.



Step 5

Unfold the Perfect Walker II so that the seat is locked into place. Push firmly down on both sides of the seat until you hear the seat click into its locked position. It will look like the image below when locked into place.



Step 6

Insert the upper arms (identified with the push pins on one side of each arm) into the right and left posts. The proper position for each is with the push pins to the outside and the thumbwheel facing forward.



Step 7

The push pins on each arm are the way to adjust the arms for an appropriate height. This can be done as a last step, if preferred.

Step 8

Insert the handles into the front of the arm tubes. You may need to loosen the thumbwheels under the arms to allow the handles to be properly inserted into the arm supports. These are adjustable. Once you have determined the appropriate length, lock the handles into place by firmly tightening the thumbwheels in the upper arms.



At this step, check all thumbwheels for tightness and that the seat is locked firmly into place. This ensures that the frame of the walker will remain open during use as well as when sitting in the seat.

The next steps are to attach the seat back onto the upper arms.

Step 9

You have been provided a set of tools along with 2 bolts, 4 washers, and 2 nuts for the attachment of the seat back.



Step 10

Identify the holes in each of the upper arms. This is where the seat back is attached to the arms.



Step 11

It is recommended that you start each side before firmly tightening. This provides some flex in the seat back and allows for an easier full install.

Step 12

Note: The curve of the seat back is positioned to the front of the walker. Also, the bolts that attach the seat strap come in from the outside of the arms so that the securing nut is on the inside of the arms (the seat side of each arm).

Line up the seat back hole to the hole in arm. Place a washer on seat strap, through the arm and then through the other side of the seat strap. Place another washer on the bolt outside of the seat strap.

Use the open ended wrench and the hex wrench to tighten the nut to the bolt. Leave it a little loose on the first side. Repeat this step for the opposite side. Now tighten both sides firmly.

Your Perfect Walker II is now ready for use. Please use it only according to the directions and instructed use guidelines.

Customer Support

If you have ANY questions, please don't hesitate to ask - you can contact our Customer Service Agents Monday - Friday between 9am and 8pm eastern time. Our toll free number is 800-958-8324.